

**VACUETTE® EasyClean**  
REF 840052

## Desinfektionsmittelliste / List of disinfection liquids

<b>Verträgliche Desinfektionstücher / compatible disinfection wipes:</b>	<b>Verträgliche Flächendesinfektionsmittel / compatible disinfection fluids:</b>
<ul style="list-style-type: none"> <li>• Ecolab Incidin™ Alcohol Wipe</li> <li>• Ecolab Incidin™ OxyWipe S</li> <li>• Hartmann BacilloI® Sensitive Tissues</li> <li>• Schülke mikroZid® universal wipes premium</li> <li>• Bode Mikrobac® Tissues</li> <li>• Dr. Schumacher Cleanisept® Wipes Maxi</li> <li>• Dr. Schumacher Ultrasol OXY® Wipes</li> <li>• B.Braun Meliseptol® sensitive Wipes</li> <li>• B.Braun Meliseptol® acute Wipes</li> <li>• Clinell® Peracetic acid wipes</li> <li>• Clinell® universal wipes</li> <li>• Sani-Cloth Disinfectant Wipe</li> </ul>	<ul style="list-style-type: none"> <li>• Bode Kohrsolin® FF</li> <li>• Ecolab Incidin™ Rapid</li> <li>• Ecolab Incidin™ Active (powder)</li> <li>• Ecolab Incidin™ Pro</li> <li>• Ecolab Sekusept™ aktiv (powder)</li> <li>• Schülke perform® (powder)</li> <li>• SoChlor® DST (1.7g NaDCC)</li> <li>• Milton™ Sterilising Tablets</li> </ul>

### Unverträgliches Flächendesinfektionsmittel / Non-compatible disinfection fluid:

- Bode BacilloI® AF