

Desinfektionsmittelliste für CBC NUCOS

List of disinfection liquids for CBC NUCOS

Verträgliche Desinfektionstücher: <i>compatible disinfection wipes:</i>	Verträgliche Flächendesinfektionsmittel: <i>compatible disinfection fluids:</i>
• Ecolab Incidin™ Alcohol Wipe • Ecolab Incidin™ OxyWipe S • Hartmann BacilloI® Sensitive Tissues • Schülke mikroZid® universal wipes premium • Bode Mikrobac® Tissues • Dr. Schumacher Cleanisept® Wipes Maxi • Dr. Schumacher Ultrasol OXY® Wipes • B. Braun Meliseptol® sensitive Wipes • B. Braun Meliseptol® acute Wipes • Clinell® Peracetic acid wipes • Clinell® universal wipes	• Bode Kohrsolin® FF • Ecolab Incidin™ Rapid • Ecolab Incidin™ Active (powder) • Ecolab Incidin™ Pro • Ecolab Sekusept™ aktiv (powder) • Schülke perform® (powder) • SoChlor® DST (1.7g NaDCC)

Unverträgliches Flächendesinfektionsmittel: <i>Non-compatible disinfection fluid:</i>
• Bode BacilloI® AF