

CBC NUCOS

The comfortable silicone tourniquet



Focus Open 2022
Special Mention

+ Comfort for Patient and Caregiver

Intuitive stasis without finger barriers or pinched skin, thanks to the silicone band with an innovative honeycomb structure. The ergonomic closure allows one-handed operation, saving caregivers valuable time – a win for comfort and efficiency.


+ Hygienic and durable

The reusable CBC NUCOS replaces single-use products, reducing waste and saving costs. Autoclavable at 134 °C, easy to disinfect with disinfection wipes, and tested for 11,200 usages and 400 sterilization cycles – a long-lasting, hygienic alternative.

+ Advantages

- + Suitable for MRI
- + Classic, sturdy snap-in closure
- + No unintentional opening in the tensioned state
- + Extra button for slow release
- + Durable, high temperature resistant material
- + Made in Germany – International patents

Learn more



+ Available in two sizes

CBC NUCOS for adults and CBC NUCOS Kids.




CBC NUCOS Kids: REF 500581
Strap: 12,5 mm x 350 mm




CBC NUCOS: REF 500580
Strap: 25 mm x 500 mm

Easy handling



1 Button for gentle release | 2 Open



1 Button for gentle release | 2 Open

How-to-use

