




CBC Penlight –
your reliable tool for
precise diagnostics.



Discover the
CBC Penlight!

Simply scan the QR code
or learn more at

www.kimetec.de



Kimetec GmbH
Dornierstr. 27
71254 Ditzingen
Germany

Tel.: +49 (0) 7156 / 17602-200
Fax.: +49 (0) 7156 / 17602-500
E-Mail: info@kimetec.de
www.kimetec.de



CBC Penlight

Diagnostic light
with LED



Efficient, durable and precise

The CBC Penlight from Kimetec ensures precise diagnostics with modern LED technology – ideal for everyday medical use and demanding applications

+ Details and advantages

- + Medical device Class I
- + Approved for eye and throat diagnostics
- + CE marking permanently laser-engraved
- + Classified according to DIN EN 62471
"Exempt Group"
- + Warm white LED – color temperature 3,000 K
- + Light intensity 15,000–31,000 mcd (millicandela)
- + Illumination time approx. 100 hours
- + LED lifetime approx. 100,000 hours
- + Cost-efficient and environmentally friendly
- + Momentary and continuous operation possible
- + Replaceable AAA batteries (incl. 2 batteries)
- + Housing made of anodized aluminum –
shatterproof and impact-resistant
- + Weight: approx. 50 g
- + Dimensions: Length 142 mm, Ø 11.5 mm
- + Made in Germany



Versatile designs for your style

The CBC Penlight is available in six colors and can be customized with your own logo as a high-quality laser engraving.



black
REF 200201



silver
REF 200250



red
REF 200270



purple
REF 200280



blue
REF 200260



lime
REF 200291

Customize
it yourself!



Also available with accessories:

Spatula Holder
REF 2002

