

mamivac® SMARLA

A new standard for **CLINICAL BREAST PUMPS**

We look forward to supporting clinics worldwide with this innovative solution.

- The mamivac® SMARLA is the answer to the requirements of clinical life
- Gentle, hygienic, silent and powerful operation guaranteed
- Intuitive touch display with digital instructions for use (IFU)
- mamivac® 2-Phase technology with both stimulation and expressing modes
- Dual program options: Alpha program is used to stimulate milk production during the first days after delivery. Pump program is used to express and collect breast milk.
- Remote control for convenient operation, especially when double-pumping
- Quiet operation for discreet use
- Flexible language options with up to 5 selectable languages
- Available with power cord or optional rechargeable battery
- Personalized settings & pump statistics stored on a chip
- Mobile version with five-footed stand optional

Made in Germany



Learn more about our
innovative products
www.kimetec.de/en

Kimetec GmbH
Gerlinger Str. 36-38
71254 Ditzingen
Germany



Available soon

mamivac® SUCTION KITS

Efficient solutions for breast milk expression for **CLINICS AND HOME USE**

mamivac® suction kits provide clinics, healthcare facilities, and mothers at home with tailored solutions for various needs. Compatible with all mamivac® breast pumps, the three BPA-free models are available in sizes S (20 mm), M (24 mm), L (28 mm), and XL (36 mm) and can be adapted to individual requirements.

mamivac® HOME with wide-neck bottle

- For home use, designed for repeated use by mothers
- Reusable for sustainable application
- Easy to use: few, easy-to-assemble parts

mamivac® CLINIC with standard bottle


Ready-to-Use: Quick, hygienic solution for clinical routines

- Germ-reduced and ready to use without prior preparation
- Suitable for up to 10 uses within 24 hours – efficient for clinical operations

Steril: For the highest hygiene standards in clinics

- Sterile packaging, immediately ready for safe use
- Single-use only – minimizes the risk of cross-contamination

Made in Germany



Learn more about our innovative products
www.kimetec.de/en

Kimetec GmbH
Gerlinger Str. 36-38
71254 Ditzingen
Germany



Available soon