Kimetec

MADE IN GERMANY


mamivac® NIPPLE SHIELDS

Protect sore and sensitive nipples

- To be used on sensitive, sore or painful nipples
- Help also with flat and retracted nipples or with baby's nipple confusion
- Ergonomic design – allows the baby to feel the mother's warmth and perceive her scent
- Best body contact between mother and child
- The well sticking nipple shields are made of thin, skin-friendly and tasteless, medical-grade silicone
- EO-sterilized or can be disinfected by boiling in water
- Medical device, Class IIa

BPA free

Conical




2 pieces with storage box

- Small (18 mm) | REF 281441
- Medium (20 mm) | REF 281448
- Large (28 mm) | REF 281442

Clinic pack (50 units loose in polybag)

- Small (18 mm) | REF 281443
- Medium (20 mm) | REF 281440
- Large (28 mm) | REF 281444

Cherry-shaped



2 pieces with storage box

- Small (18 mm) | REF 281481
- Medium (22 mm) | REF 281482
- Large (28 mm) | REF 281515

Clinic pack (50 units loose in polybag)

- Small (18 mm) | REF 281483
- Medium (22 mm) | REF 281484
- Large (28 mm) | REF 281516

Conical (sterile)

Individually blister-packed
(50 pieces loose in carton)

- Small (18 mm)
  REF 281601
- Medium (20 mm)
  REF 281604
- Large (28 mm)
  REF 281607

Storage box

25 pieces in box

- REF 281446



CE 0124

Kimetec

MADE IN GERMANY


mamivac® MAMILLOMETER

To find the correct size of nipple and breast shields

To determine the correct nipple shield size, we recommend using our Mamillometer. It is essential that the size of the nipple shield is based on the mamilla (nipple) rather than the areola. The correct size is achieved when there is approximately 2 mm of space around the nipple.

Sizing pad with 15 sheets

- REF WT1000057GB



mamivac® DEMO BOX

All variants in practical folding box for demonstration

Box with 6 pieces

- REF 281447



Learn more about our innovative products
www.kimetec.de/en

#welovemedtech

Kimetec GmbH
Dornierstr. 27
71254 Ditzingen
Germany

P281448EN_2026-04